



Lab Project Enrichment & N Cycle

Enrichment Techniques

Goal - tease target microbes from population
Inoculum - mixed sample with target microbes
Many strategies thereafter - physical, chemical
Understand inoculum, strategies - ALL projects!

Nitrogen Fixation

N₂ fixed to NH₃ and NH₂ groups - COSTS ATP
ONLY by Prokaryotes - Bacteria and Archaea
Conserved O₂-intolerant Nitrogenase/Fe-Mo
Understand Acetylene Reduction Assay

Some Nitrogen Fixer Diversity

Cyanobacteria: colonial, heterocysts resist O₂
Methanogens/Archaea: live in anoxic habitats
PJ - Rhizobium: symbiotic leghemoglobin
PJ - Azotobacter: fastest ETC uses O₂

Symbiotic Nitrogen Fixation

Legumes: root hairs secrete inducer molecules
Symbionts: soil Rhizobium, Proteobacteria
Separate/no fixation; together/leghemoglobin
Bacteroids: misshapen Rhizobium in nodules

Free-Living Aerobic Azotobacter

Proteobacteria - capsule-producing rods
Beijerinck - first historic enrichment (Fig. 22.1)

Nitrification, Aerobic Lithotrophy

Two-step - ammonia then nitrite oxidation
NH₃ + O₂ yields NO₂ + H₂O, Nitrosomonas
NO₂ + O₂ yields NO₃, Nitrobacter
Both are Proteobacteria; packed with lamellae

Anaerobic Nitrate Reduction

Denitrifying (to N₂) OR Nitrifying (to N-fixed)
NO₃ + CH₂O yields N₂ + CO₂, Pseudomonas
NO₃ + CH₂O yields NO₂, NH₃ + CO₂, Enterics

ACTIVITIES

Azotobacter - Work Individually, Incubate RT°

Two N-free enrichment flasks have been selecting for 3 weeks; everyone needs to share!
Looping from oxic TOP of enrichment, streak onto 1 N-free plate - watch for slimy colonies

Notebook records – thoroughly describe plates, target colony texture and shape

Azotobacter N-Free Media

<i>Solution A</i>	<i>Solution B</i>
K ₂ HPO ₄ : 1.6 g	MgSO ₄ : 0.4 g
KH ₂ PO ₄ : 0.4 g	CaSO ₄ : 0.2 g
	FeSO ₄ /7H ₂ O: 0.006 g
	MoO ₃ : 0.002 g
	sucrose: 10 g

Prepared as liquid and agar (15-20 g/L) plates, with agar added to Solution B. Combine 1A:1B after autoclaving. This high-sugar enrichment medium is designed to enhance capsule production, a distinct Azotobacter trait. This medium is also highly selective for N₂-fixers because it contains no nitrogen.

Nitrate Reduction - Work In Pairs, Incubate RT°

Two 10-ml, 10-fold organic soil dilutions have been prepared; everyone needs to share!
Working in pairs, inoculate 3 nitrate reduction tubes: soil dilution, E. coli, Pseudomonas
After a week, complete all portions of the nitrate reduction test as described below

Notebook records – thoroughly describe tubes and interpret testing at each stage

Nitrate Reduction Media

Beef Extract: 3 g
Peptone: 5 g
KNO₃: 1 g

Prepared as 7 ml light gold liquid tubes with Durham tubes. After growth, read in order: (1) Bubble? (yes = NO₃ to N₂); (2) Add 10 drops A & B. Red? (yes = NO₂); (3) If not, add 10 grains Zn and wait 2-3 minutes. Red? Yes = NO₃ present/no NO₃ reduction; No = NH₃, NH₂-compounds.

Rhizobium - Work In Pairs, Incubate RT°

Working in pairs, obtain 2 root nodules, lightly wash off dirt and surface treat 2 minutes, comparing...
1 nodule with bleach disinfectant; 1 nodule with alcohol antiseptic – KEEP TRACK which is which
After treatment, separately place each nodule in 1 drop sterile water - crush with sterilized forceps
PROPERLY streak each treatment onto provided Rhizobium plates – watch for slimy colonies!

Notebook records - thoroughly describe and compare plates/treatments, target colonies

Rhizobium Media

Mannitol : 10 g
Yeast Extract: 1.0 g
MgSO₄/7H₂O: 0.2 g
NaCl: 0.2 g
K₂HPO₄: 0.5 g
FeCl₃: 0.005 g

Prepared as agar plates (15-20 g/L). This high-sugar medium is designed to enhance capsule production, a distinctive Rhizobium trait. These plates are not particularly selective and must be used in conjunction with surface-sterilization for enrichment purposes.

LAB MATERIALS TURN-IN

11 pts. Informal Notebook records: make sure all guidelines in syllabus and above are followed.

2 pts. INDIVIDUAL TURN-IN: Each person turns in his/her final Azotobacter plate, additionally CIRCLING an Azotobacter colony. It will be graded on labeling, pattern, Azotobacter isolation/ID.

2 pts. PAIR TURN-IN: Each pair turns in their set of Rhizobium plates, circling a Rhizobium colony. It will be graded as above; make sure I can tell how each was treated!